

Lecture 2

Visualization of cells

Outline

1. Light microscope
2. Fluorescence microscope
3. Confocal microscope
4. nanoscope
5. FRAP/FRET/Photoactivation
6. Electron microscope

Microscopes




Figure 9-2 Molecular Biology of the Cell 6e (© Garland Science 2015)

1. Light microscope

Bright field



(A)

Phase contrast




(B)



(C)

Differential-interference-contrast
(DIC)




(D)

dark field

50 μm

Conventional Light microscope has a resolution of $0.2\mu\text{m}$

LENSES



RESOLUTION: the resolving power of the microscope depends on the width of the cone of illumination and therefore on both the condenser and the objective lens. It is calculated using the formula

$$\text{resolution} = \frac{0.61 \lambda}{n \sin \theta}$$

where:

θ = half the angular width of the cone of rays collected by the objective lens from a typical point in the specimen (since the maximum width is 180° , $\sin \theta$ has a maximum value of 1)

n = the refractive index of the medium (usually air or oil) separating the specimen from the objective and condenser lenses


λ = the wavelength of light used (for white light a figure of $0.53 \mu\text{m}$ is commonly assumed)

The smaller the number is, the bigger the resolution is

NUMERICAL APERTURE: $n \sin \theta$ in the equation above is called the numerical aperture of the lens (NA) and is a function of its light-collecting ability. For dry lenses this cannot be more than 1, but for oil-immersion lenses it can be as high as 1.4. The higher the numerical


aperture, the greater the resolution and the brighter the image (brightness is important in fluorescence microscopy). However, this advantage is obtained at the expense of very short working distances and a very small depth of field.

1.1. Bright field microscope




Magnification : up to 1000 fold
Eye piece: up to 10X
Objective: up to 100X

1.2. Dark field microscopy



Light illuminates on the surface of the specimen



An opaque blocker to block lights directly penetrating specimen

Dark field is especially useful for finding cells in suspension

Comparison for Dark-field and Bright-field microscopes


- Dark field: cells bright, background dark.

Light illuminates from the side and only scattered light enter lenses

- Bright field: cells dark, background bright.

light pass through a cell to form image directly


1.3. Phase contrast



Phase contrast estimates phase changes after light passes through parts of a cell with different density, and then convert the change of phase into different brightness. Phase contrast microscope is frequently used in Cell Biology


Phase difference results in difference in brightness

TWO WAVES IN PHASE




BRIGHT

TWO WAVES OUT OF PHASE



DIM


Structure of a Phase Contrast Microscope



1.4. Differential-Interference Contrast (DIC)

- Analyze the optical path length **gradients** (the rate of change for refractive light).
- No annular diaphragm and phase plate needed
- Produce shadow-cast image.
- Elimination of halo artifacts
- Excellent images with relatively thick specimens


2. Fluorescence microscope





(A)

- ◆ A fluorophore absorbs a photon to be excited
- ◆ It is not stable and undergoes conformational change to emit a photon
- ◆ After it emits, it goes back to the ground state.
- ◆ The cycle can be repeated.
- ◆ The energy of excitation light is higher than the energy of emission light.
- ◆ Lower energy = longer wavelength
- ◆ Stokes shift: the difference between excitation energy and the emission energy for a given fluorophore

How does the fluorescence microscopy work?



2.1. Fluorescence probes




2.2. Categories for Fluorescence microscope

- Direct immunofluorescence primary antibody-conjugated fluorescence
- Indirect immunofluorescence secondary antibody-conjugated fluorescence




primary antibody:
rabbit antibody
directed against
antigen A


immobilized
antigen A



secondary antibodies:
marker-coupled antibodies
directed against rabbit
antibodies



A Multi-fluorescent cell



10 μm

- Step 1: cell fixation (2-4% paraformaldehyde,
acetone/methanol(1:1)
10% formalin/10%methanol)
- Step 2: permeabilize with detergent
(0.3%triton-X-100)
- Step 3: incubate with primary antibody
- Step 4: incubate with secondary-Ab-
conjugated-fluorescence

1. Can be protein-specific
2. Detect cellular localization
3. Roughly quantify a biomolecule

Aequorin, a luminescent protein that glows blue in the presence of Ca^{2+}






Figure 9-30 Molecular Biology of the Cell 6e (© Garland Science 2015)

Fura-2 can sense free Ca^{2+} ions in cells




Red: higher Ca^{2+}
Blue: lower Ca^{2+}

3. Confocal microscope



1. More accurately detect protein localization in a 3-D location.
2. Use laser beam as light source provide better resolution.
3. Use focused light from a pinhole to illuminate on a certain spot in specimen.
4. Out of focus light will be excluded

Conventional immunofluorescence



(A)


confocal microscope



(B)

10 μm

3.1. Laser scanning confocal




one pinhole

Good: better resolution

Bad: Illumination on specimen
for longer time, sometimes
cause bleaching for fluorescence signal

3.2. Spinning disk confocal



Multiple pinholes allow specimen
To be illuminated at multiple
Points simultaneously

Good: quick and dynamic imaging
Bad: lower resolution

4 Nanoscopy

Fluorescence microscopy that surpasses the limitation of $0.2\mu\text{m}$



Stefan W. Hell

Prize share: 1/3



William E. Moerner



Prize share: 1/3




Eric Betzig

Prize share: 1/3

Nobel Laureates
in 2014 Chemistry



Caveats with Super-Resolution Microscopy: Live Cells




Nyquist criterion:

N -fold resolution increase in D dimensions $\rightarrow N^D$ -fold more photons collected

L. Schermelleh, R. Heintzmann, *J. Cell Biol.* (2010)

STED / RESOLFT



reported
resolution
(nm)

xy: 20 nm

xyz: 30 nm

photon
increase
required

100

1,070

intensity
(W/cm²)


$10^4 - 10^9$

acquisition
time (sec)

> 60

~1,000

Localization



xy: 20 nm

xy: 10 nm,
z: 20 nm

100


14,400

$10^3 - 10^4$

>20

1,500

SIM



xy: 100 nm

xy: 100 nm,
z: 370 nm

4




8

$10 - 10^2$

0.1 - 1



~10

The Beauty and Complexity of Living Systems



--- Nobel prize website

5.1 FRET (fluorescence resonance energy transfer) to detect protein-protein interaction .




Two separate dyes to label two interesting molecules, respectively.

Upon direct interaction, The energy from one dye will cause the excitation of the other dye and to emit light. The signal can be Recorded.



(B) NO PROTEIN INTERACTION



NO EXCITATION OF GREEN
FLUORESCENT PROTEIN,
BLUE LIGHT DETECTED




(C) PROTEIN INTERACTION

FLUORESCENCE RESONANCE
ENERGY TRANSFER,
GREEN LIGHT DETECTED

e.g. FRET



(D)




(E)

FRET occurs in close distance
~5nm.



Rules for two fluorochromes:
Emission spectrum of one
should overlap with the
absorption spectrum of the
other.

5.2 Photoactivation



Activate a target molecule by photoactivation, as time goes, the fluorescence signal changes. It can measure protein trafficking, turnover, degradation, diffusion, etc.


5.3 Fluorescence recovery after photobleaching (FRAP)





10.6-FRAP.mov

6. Electron microscope




Resolution: 0.1nm-1nm


Transmission microscope

Special sample preparation

- Sample very thin: 50-100nm thick(penetrating power for electron is very weak)
- Sample rapidly frozen to preserve its original state
- Glutaraldehyde(to crosslink proteins) and osmium tetroxide (to stabilize lipid membrane)
- Sample dehydrated, resin-embedded




glutaraldehyde



osmium tetroxide

An image of a cell under electron microscope




Electron microscope

- Immunogold electron microscope
- Scanning electron microscope (SEM)
- Transmission electron microscope (TEM)

6.1. Immunogold electron microscope

- Special protein is incubated with primary antibody
- Gold particle attached secondary antibody is incubated to label the protein with a gold particle.
- Gold is electron-dense and can be seen under the electron microscope


Images of immunogold electron microscope



6.2. Scanning electron microscope

- Produce 3-D structure for the surface of a specimen.
- Uses electrons that are scattered or emitted from the specimen's surface.
- Lower resolution ~ 10nm
- Usually whole cells or tissues, not subcellular organelles.
- Less expensive


SEM



6.3. Transmission electron microscope

- Electrons penetrate through the specimen
- Produce high resolution images
- Observed subcellular organelles, virus structure, etc.
- More expensive

TEM



(b) Transmission electron microscope

0.5 micrometers